Atty. Docket No. BI-15 Combine Mobile Phase Flow Stream in Chromatography System Separate Liquid Phase Fractions and Send to Collection System Provide Extended Vessel Assembly for Collection System Receive Liquid Phase Through Mouth of Vessel Extender Fill Assembly With Liquid Phase Beyond Volume of Collection Vessel Provide Volume in the Vessel Extender for Expansion of Gas Phase Without Disturbing the Surface of Liquid Phase to Point of Aerosol Formation Reduce Liquid Volume Using **Dry Down Process**

NEWHEPLACEMENT SHEET

Application No.: 10/705,907




FIGURE 2